

## PRK2 antibody

Catalog No: #22595

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PRK2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 46 and 336 of Human PKN2
Target Name	PRK2
Accession No.	Swiss-Prot:Q16513Gene ID:5586
Uniprot	Q16513
GeneID	5586;
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 112kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

## Images




Sample (30 ug of whole cell lysate)


A: H1299

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using protein kinase N2 antibody at 1: 200 dilution.

---

Note: This product is for in vitro research use only