

## Phenethyl alcohol


Catalog No: #S0011

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                   |   |
|-------------------|---|
| Product Name      | Phenethyl alcohol                               |
| Brief Description | Inhibitors                                      |
| Purification      | 99.82%  |
| Target Name       | Others  |
| Calculated MW     | 122.16  |
| Formulation       | C <sub>8</sub> H <sub>10</sub> O                |
| Storage           | 3 years -20°C powder; 2 years -80°C in solvent; |

## Images



## Product Description

Research Area: Immunology SMILES: OCCC1=CC=CC=C1 Pathways: Others Receptor: Others Boiling pt: 219-221°C Melting pt: -27°C Solubility: soluble in Water Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Lingappa BT, et al. Science. 1969 Jan 10;163(3863):192-4.

Note: This product is for in vitro research use only