

CLC-4 Polyclonal Antibody

Catalog No: #40755

Package Size: #40755-1 50ul #40755-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CLC-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	CLC-4 Polyclonal Antibody detects endogenous levels of CLC-4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CLC-4.
Target Name	CLC-4
Other Names	CLCN4; H(+)/Cl(-) exchange transporter 4; Chloride channel protein 4; CIC-4; Chloride transporter CIC-4
Accession No.	Swiss-Prot: P51793NCBI Gene ID: 1183
Uniprot	P51793
GeneID	1183;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of RAT-MUSCLE cells using CLC-4 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using CLC-4 Polyclonal Antibody

Note: This product is for in vitro research use only