

CLC-6 Polyclonal Antibody

Catalog No: #40756

Package Size: #40756-1 50ul #40756-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLC-6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt Mk
Specificity	CLC-6 Polyclonal Antibody detects endogenous levels of CLC-6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CLC-6.
Target Name	CLC-6
Other Names	CLCN6; KIAA0046; Chloride transport protein 6; Chloride channel protein 6; CIC-6
Accession No.	Swiss-Prot: P51797NCBI Gene ID: 1185
Uniprot	P51797
GeneID	1185;
SDS-PAGE MW	97kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of L929 cells using CLC-6 Polyclonal Antibody



Western Blot analysis of K562 cells using CLC-6 Polyclonal Antibody

Note: This product is for in vitro research use only