

Cephadrine


Catalog No: #S0199

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cephadrine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibacterial inhibitor
Calculated MW	349.41
Formulation	C16H19N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=C(N2[C@@H]([C@@H](C2=O)NC(=O)[C@@H](C3=CCC=CC3)N)SC1)C(=O)O Pathways: Microbiology & virology Receptor: PBPs Boiling pt: 693.1°C Melting pt: 140-142°C Solubility: soluble in Water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sharma UK, et al. J Bacteriol. 1995 Dec; 177(23):6745-50.

Note: This product is for in vitro research use only