

Chromocarbe


Catalog No: #S0259

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chromocarbe
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	190.15
Formulation	C10H6O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)C1=CC(=O)C2=CC=CC=C2O1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 260°C
Solubility: DMSO: 38 mg/mL (199.84 mM)
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Godeau G, et al. Clin Physiol Biochem, 1987, 5(1), 15-26.

Note: This product is for in vitro research use only