

Isoxicam


Catalog No: #S0285

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isoxicam
Brief Description	Inhibitors
Purification	98.00%
Target Name	COX inhibitor
Calculated MW	335.34
Formulation	C ₁₄ H ₁₃ N ₃ O ₅ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CN1C(C(=O)NC2=NOC(C)=C2)=C(O)C2=CC=CC=C2S1(=O)=O
Pathways: Immunology/Inflammation;
Neuroscience
Receptor: COX
Boiling pt: 518.8°C
Melting pt: 265-271°C
Solubility: DMSO: 10 mM
Appearance: Faint Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Burch FX. Am J Med. 1985 Oct 18;79(4B):28-32.

Note: This product is for in vitro research use only