

COL11A1 Polyclonal Antibody

Catalog No: #40770



Package Size: #40770-1 50ul #40770-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	COL11A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	COL11A1 Polyclonal Antibody detects endogenous levels of COL11A1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human COL11A1.
Target Name	COL11A1
Other Names	COL11A1; COLL6; Collagen alpha-1(XI) chain
Accession No.	Swiss-Prot: P12107NCBI Gene ID: 1301
Uniprot	P12107
GeneID	1301;
SDS-PAGE MW	181kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of KB cells using COL11A1 Polyclonal Antibody

Note: This product is for in vitro research use only