

Xylocaine


Catalog No: #S0468

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Xylocaine
Brief Description	Inhibitors
Purification	99.74%
Target Name	Histamine Receptor agonist
Calculated MW	234.34
Formulation	C14H22N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: CCN(CC)CC(=O)NC1=CC(C)=CC=C1C
Pathways: GPCR/G Protein; Immunology/Inflammation; Neuroscience
Receptor: H1 receptor
Boiling pt: Melting pt: 66-69°C
Solubility: Ethanol: 47 mg/mL (200.56 mM); Water: 9 mg/mL (38.4 mM); DMSO: 47 mg/mL (200.56 mM)
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sheets PL, et al. J Physiol. 2007 Jun 15;581(Pt 3):1019-31.

Note: This product is for in vitro research use only