

L-Pyroglutamic acid


Catalog No: #S0507

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-Pyroglutamic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	129.11
Formulation	C5H7NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:C1CC(=O)N[C@@H]1C(=O)OPathways:OthersReceptor:Unusual Amino AcidsBoiling pt:Melting pt:160-163CSolubility:soluble in WaterAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Imming P, et al. Nat Rev Drug Discov. 2006 Oct; 5(10):821-34.

Note: This product is for in vitro research use only