

3-Methyloxindole


Catalog No: #S0538

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-Methyloxindole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	147.18
Formulation	C ₉ H ₉ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC1C(=O)NC2=CC=CC=C12
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Potchoiba MJ, et al. Am J Vet Res. 1982 Aug; 43(8):1418-23.

Note: This product is for in vitro research use only