

Hydrocaffeic acid


Catalog No: #S0539

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hydrocaffeic acid
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	182.18
Formulation	C ₉ H ₁₀ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)CCC1=CC(O)=C(O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Off-white Color Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Aoki R, et al. J Neurosci. 2011 Nov 16; 31(46): 16603-10.

Note: This product is for in vitro research use only