

Trimethylacetic anhydride


Catalog No: #S0546

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Trimethylacetic anhydride
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	186.25
Formulation	C ₁₀ H ₁₈ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC(C)(C)C(=O)OC(=O)C(C)(C)C
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Shiina I, et al. Chemistry. 2010 Jan 4;16(1):167-72.

Note: This product is for in vitro research use only