

Domperidone EP Impurity A

Catalog No: #S0550

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Domperidone EP Impurity A
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	251.71
Formulation	C ₁₂ H ₁₄ ClN ₃ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: ClC1=CC2=C(C=C1)N(C1CCNCC1)C(=O)N2 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only