

C-TAK1 Polyclonal Antibody

Catalog No: #40793



Package Size: #40793-1 50ul #40793-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-TAK1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	C-TAK1 Polyclonal Antibody detects endogenous levels of C-TAK1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human C-TAK1.
Target Name	C-TAK1
Other Names	MARK3; CTAK1; EMK2; MAP/microtubule affinity-regulating kinase 3; C-TAK1; cTAK1; Cdc25C-associated protein kinase 1; ELKL motif kinase 2; EMK-2; Protein kinase STK10; Ser/Thr protein kinase PAR-1; Par-1a; Serine/threonine-protein kinase p78
Accession No.	Swiss-Prot: P27448NCBI Gene ID: 4140
Uniprot	P27448
GeneID	4140;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using C-TAK1 Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN cells using C-TAK1 Polyclonal Antibody

Note: This product is for in vitro research use only