

Diphenyl sulfide


Catalog No: #S0574

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diphenyl sulfide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	186.27
Formulation	C ₁₂ H ₁₀ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: S(C1=CC=CC=C1)C1=CC=CC=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Michida T, et al. Yakugaku Zasshi. 1999 Oct;119(10):780-5.

Note: This product is for in vitro research use only