

CUL-2 Polyclonal Antibody

Catalog No: #40797

Package Size: #40797-1 50ul #40797-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CUL-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	CUL-2 Polyclonal Antibody detects endogenous levels of CUL-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human CUL-2.
Target Name	CUL-2
Other Names	CUL2; Cullin-2; CUL-2
Accession No.	Swiss-Prot: Q13617NCBI Gene ID: 8453
Uniprot	Q13617
GeneID	8453;
SDS-PAGE MW	86kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using CUL-2 Polyclonal Antibody



Western Blot analysis of Jurkat NIH-3T3 K562 cells using CUL-2 Polyclonal Antibody

Note: This product is for in vitro research use only