

Sulisobenzone


Catalog No: #S0691

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sulisobenzone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	308.31
Formulation	C ₁₄ H ₁₂ O ₆ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=C(C=C(C(=O)C2=CC=CC=C2)C(O)=C1)S(=O)(=O)O
Pathways: Others
Receptor: Others
Boiling pt: 0°C
Melting pt: 170°C
Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nohynek GJ, Schaefer H. Regul Toxicol Pharmacol. 2001 Jun;33(3):285-99.

Note: This product is for in vitro research use only