

CYP26A1 Polyclonal Antibody

Catalog No: #40810

Package Size: #40810-1 50ul #40810-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CYP26A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	CYP26A1 Polyclonal Antibody detects endogenous levels of CYP26A1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CYP26A1.
Target Name	CYP26A1
Other Names	CYP26A1; CYP26; P450RAI1; Cytochrome P450 26A1; Cytochrome P450 retinoic acid-inactivating 1; Cytochrome P450RAI; hP450RAI; Retinoic acid 4-hydroxylase; Retinoic acid-metabolizing cytochrome
Accession No.	Swiss-Prot: O43174NCBI Gene ID: 1592
Uniprot	O43174
GeneID	1592;
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of 293 cells using CYP26A1 Polyclonal Antibody



Western Blot analysis of HepG2 cells using CYP26A1 Polyclonal Antibody

Note: This product is for in vitro research use only