

CYP2S1 Polyclonal Antibody

Catalog No: #40816

Package Size: #40816-1 50ul #40816-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CYP2S1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	CYP2S1 Polyclonal Antibody detects endogenous levels of CYP2S1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CYP2S1.
Target Name	CYP2S1
Other Names	CYP2S1; Cytochrome P450 2S1; CYP11S1
Accession No.	Swiss-Prot: Q96SQ9NCBI Gene ID: 29785
Uniprot	Q96SQ9
GeneID	29785;
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using CYP2S1 Polyclonal Antibody

Note: This product is for in vitro research use only