

Butyl 4-aminobenzoate


Catalog No: #S0916

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Butyl 4-aminobenzoate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Sodium Channel inhibitor
Calculated MW	193.24
Formulation	C ₁₁ H ₁₅ NO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCCCOC(=O)C1=CC=C(N)C=C1
Pathways: Membrane transporter/Ion channel
Receptor: Sodium Channel
Boiling pt: 174°C
Melting pt: 57-58°C
Solubility: DMSO: 10 mM
Appearance: White Needle Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Tzimopoulos D, et al. J Inorg Biochem. 2010 Apr;104(4):423-30.

Note: This product is for in vitro research use only