

Docusate sodium

Catalog No: #S0929


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Docusate sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Opioid Receptor antagonist
Calculated MW	444.56
Formulation	C ₂₀ H ₃₇ NaO ₇ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: [Na+].CCCCC(CC)COC(=O)CC(C(=O)OCC(CC)CCCC)S([O-])(=O)=O Pathways: Endocrinology/Hormones; GPCR/G Protein;

Neuroscience Receptor: Opioid Receptor Boiling pt: Melting pt: 173-179°C Solubility: Soluble in Water Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nader ME, et al. Otolaryngol Head Neck Surg. 2012 Mar;146(3):455-60.

Note: This product is for in vitro research use only