

DDX8 Polyclonal Antibody

Catalog No: #40838

Package Size: #40838-1 50ul #40838-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DDX8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	DDX8 Polyclonal Antibody detects endogenous levels of DDX8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human DDX8.
Target Name	DDX8
Other Names	DHX8; DDX8; ATP-dependent RNA helicase DHX8; DEAH box protein 8; RNA helicase HRH1
Accession No.	Swiss-Prot: Q14562NCBI Gene ID: 1659
Uniprot	Q14562
GeneID	1659;
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using DDX8 Polyclonal Antibody

Note: This product is for in vitro research use only