

Dynactin 2 Polyclonal Antibody

Catalog No: #40856

Package Size: #40856-1 50ul #40856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dynactin 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Dynactin 2 Polyclonal Antibody detects endogenous levels of Dynactin 2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Dynactin 2.
Target Name	Dynactin 2
Other Names	DCTN2; DCTN50; Dynactin subunit 2; 50 kDa dynein-associated polypeptide; Dynactin complex 50 kDa subunit; DCTN-50; p50 dynamitin
Accession No.	Swiss-Prot: Q13561NCBI Gene ID: 10540
Uniprot	Q13561
GeneID	10540;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using Dynactin 2 Polyclonal Antibody

Note: This product is for in vitro research use only