

# Dynein LC 1 Polyclonal Antibody

Catalog No: #40857



Package Size: #40857-1 50ul #40857-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | Dynein LC 1 Polyclonal Antibody                                                                                       |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB IHC ELISA                                                                                                          |
| Species Reactivity    | Hu Ms                                                                                                                 |
| Specificity           | Dynein LC 1 Polyclonal Antibody detects endogenous levels of Dynein LC 1 protein.                                     |
| Immunogen Type        | peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from the Internal region of human Dynein LC 1.                                            |
| Target Name           | Dynein LC 1                                                                                                           |
| Other Names           | DNAL1; C14orf168; Dynein light chain 1; axonemal                                                                      |
| Accession No.         | Swiss-Prot: Q4LDG9NCBI Gene ID: 83544                                                                                 |
| Uniprot               | Q4LDG9                                                                                                                |
| GeneID                | 83544;                                                                                                                |
| SDS-PAGE MW           | 22kd                                                                                                                  |
| Concentration         | 1mg/ml                                                                                                                |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of MCF7 cells using Dynein LC 1 Polyclonal Antibody

---

Note: This product is for in vitro research use only