

Dyrk1A Polyclonal Antibody

Catalog No: #40858

Package Size: #40858-1 50ul #40858-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dyrk1A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Dyrk1A Polyclonal Antibody detects endogenous levels of Dyrk1A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Dyrk1A.
Target Name	Dyrk1A
Other Names	DYRK1A; DYRK; MNB; MNBH; Dual specificity tyrosine-phosphorylation-regulated kinase 1A; Dual specificity YAK1-related kinase; HP86; Protein kinase minibrain homolog; MNBH; hMNB
Accession No.	Swiss-Prot: Q13627NCBI Gene ID: 1859
Uniprot	Q13627
GeneID	1859;
SDS-PAGE MW	90kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Dyrk1A Polyclonal Antibody



Western Blot analysis of VEC cells using Dyrk1A Polyclonal Antibody

Note: This product is for in vitro research use only