

Oxiconazole nitrate

Catalog No: #S1290


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Oxiconazole nitrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Antifection; Others
Calculated MW	492.14
Formulation	C18H13Cl4N3OHNO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1=CC(=C(C=C1Cl)Cl)CO/N=C(\CN2C=CN=C2)/C3=C(C=C(C=C3)Cl)Cl.[N+](=O)(O)[O-] Pathways: Microbiology&virology;
 Others Receptor: 14-α Demethylase; LSS Boling pt: 576.8C Melting pt: 137-138C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Carrillo-Munoz AJ, et al. Rev Esp Quimioter. 2006 Jun;19(2):130-9.

Note: This product is for in vitro research use only