

Tiapride hydrochloride


Catalog No: #S1292

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tiapride hydrochloride
Brief Description	Inhibitors
Purification	99.64%
Target Name	Dopamine Receptor antagonist
Calculated MW	364.89
Formulation	C15H24N2O4SHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: Cl.CCN(CC)CCNC(=O)C1=C(OC)C=CC(=C1)S(C)(=O)=O
Pathways: GPCR/G Protein;
Neuroscience
Receptor: D2; D3
Boling pt: 498.1°C
Melting pt: Solubility: Water: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Scatton B, et al. Eur Psychiatry. 2001 Jan;16 Suppl 1:29s-34s.

Note: This product is for in vitro research use only