

# Sucralfate


Catalog No: #S1338

Orders: order@signalwayantibody.com  
Support: tech@signalwayantibody.com

## Description

Product Name	Sucralfate
Brief Description	Inhibitors
Purification	98.00%
Target Name	FGFR; Others
Calculated MW	2086.75
Formulation	C12H54Al16O75S8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism

systemSMILES:O.[Al].[Al].[Al].[Al].[Al].[Al].[Al].[Al].[Al].OS(=O)(=O)OC[C@H]1O[C@@H](O[C@]2(COS(O)(=O)O[C@H](OS(O)(=O)=O)[C@@H](OS(O)(=O)=O)[C@@H]2OS(O)(=O)=O)[C@H](OS(O)(=O)=O)[C@@H](OS(O)(=O)=O)[C@@H]1OS(O)(=O)=O)OP

Pathways:Angiogenesis; Others; Tyrosine Kinase/AdaptorsReceptor:FGF2; Fibrinogen chain; Pepsin; PGFBoling pt:Melting pt:Solubility:Insoluble in cold Water.Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Terao N, et al. Arzneimittelforschung. 1980;30(1):76-8.

Note: This product is for in vitro research use only