

Methyldienedione


Catalog No: #S1377

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methyldienedione
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor antagonist
Calculated MW	270.37
Formulation	C18H22O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CCC3=C4CCC(=O)C=C4CC[C@H]3[C@@H]1CCC2=OPathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Scarth JP, et al. Steroids. 2010 Oct;75(10):643-52.

Note: This product is for in vitro research use only