

EKLF/CKLF/UKLF Polyclonal Antibody

Catalog No: #40873



Package Size: #40873-1 50ul #40873-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EKLF/CKLF/UKLF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	EKLF/CKLF/UKLF Polyclonal Antibody detects endogenous levels of EKLF/CKLF/UKLF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human EKLF/CKLF/UKLF.
Target Name	EKLF/CKLF/UKLF
Other Names	KLF1; EKLF; Krueppel-like factor 1; Erythroid krueppel-like transcription factor; EKLF; KLF5; BTEB2; CKLF; IKLF; Krueppel-like factor 5; Basic transcription element-binding protein 2; BTE-binding protein 2; Colon krueppel-like factor; GC-bo
Accession No.	Swiss-Prot: Q13351/Q13887/O75840NCBI Gene ID: 10661/688/8609
Uniprot	Q13351
GeneID	10661;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of VEC cells using EKLF/CKLF/UKLF Polyclonal Antibody



Western Blot analysis of VEC HY926 cells using EKLF/CKLF/UKLF Polyclonal Antibody

Note: This product is for in vitro research use only