

Elf-1 Polyclonal Antibody

Catalog No: #40874

Package Size: #40874-1 50ul #40874-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elf-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Elf-1 Polyclonal Antibody detects endogenous levels of Elf-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Elf-1.
Target Name	Elf-1
Other Names	ELF1; ETS-related transcription factor Elf-1; E74-like factor 1
Accession No.	Swiss-Prot: P32519NCBI Gene ID: 1997
Uniprot	P32519
GeneID	1997;
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Elf-1 Polyclonal Antibody

Note: This product is for in vitro research use only