

NCKAP1 antibody

Catalog No: #22609

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	NCKAP1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 699 and 912 of Human NCKAP1
Target Name	NCKAP1
Accession No.	Swiss-Prot:Q9Y2A7Gene ID:10787
Uniprot	Q9Y2A7
GeneID	10787;
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 130kd

Western blotting: 1:500-1:3000

Images



Sample (30ug whole cell lysate)
A: Hep G2
7.5% SDS PAGE
Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only