

EphB4 Polyclonal Antibody

Catalog No: #40891

Package Size: #40891-1 50ul #40891-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EphB4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	EphB4 Polyclonal Antibody detects endogenous levels of EphB4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human EphB4.
Target Name	EphB4
Other Names	EPHB4; HTK; MYK1; TYRO11; Ephrin type-B receptor 4; Hepatoma transmembrane kinase; Tyrosine-protein kinase TYRO11
Accession No.	Swiss-Prot: P54760NCBI Gene ID: 2050
Uniprot	P54760
GeneID	2050;
SDS-PAGE MW	108kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HY926 cells using EphB4 Polyclonal Antibody

Note: This product is for in vitro research use only