

Chlormadinone acetate


Catalog No: #S1596

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chlormadinone acetate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Androgen Receptor antagonist
Calculated MW	404.93
Formulation	C23H29ClO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Endocrine

system SMILES: CC(=O)[C@]1(CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2C=C(C4=CC(=O)CC[C@]34C)Cl)C)OC(=O)C Pathways: Endocrinology/Hormones
Receptor: Androgen Receptor
Boiling pt: 512.5°C
Melting pt: 212°C
Solubility: Water: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kim JM, et al. Eur J Pharmacol. 2014 Mar 5;726:1-8.

Note: This product is for in vitro research use only