

Actriol


Catalog No: #S1683

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Actriol
Brief Description	Inhibitors
Purification	98.00%
Target Name	GR inhibitor
Calculated MW	288.39
Formulation	C18H24O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Endocrine

system SMILES: C[C@]12CC[C@H]3[C@H]([C@@H]1C[C@@H]([C@@H]2O)CCC4=C3C=CC(=C4)O Pathways: Metabolism Receptor: Glucocorticoid Receptor
Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Miller E, et al. J Interferon Cytokine Res. 1998 Nov;18(11):921-5.

Note: This product is for in vitro research use only