


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TG003
Brief Description	Inhibitors
Purification	99.04%
Target Name	BCL inhibitor
Calculated MW	249.33
Formulation	C13H15NO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCN1C2=C(C=CC(=C2)OC)S/C1=C\C(=O)C
Pathways: Apoptosis
Receptor: mCLK1; mCLK2; mCLK4
Boiling pt: Melting pt: Solubility: DMSO: 6 mg/mL warmed (24.06 mM)
Appearance: Off-white Color Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Muraki M, et al. J Biol Chem. 2004, 279(23), 24246-24254.

Note: This product is for in vitro research use only