

Dabrafenib (GSK2118436)


Catalog No: #S1903

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dabrafenib (GSK2118436)
Brief Description	Inhibitors
Purification	99.50%
Target Name	Raf inhibitor
Calculated MW	519.56
Formulation	C23H20F3N5O2S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: CC(C)(C)C1=NC(=C(S1)C1=NC(N)=NC=C1)C1=C(F)C(NS(=O)(=O)C2=C(F)C=CC=C2F)=CC=C1 Pathways: MAPK Receptor: B-Raf; B-Raf (V600E); C-Raf
Boiling pt: Melting pt: Solubility: DMSO: 30 mg/mL (57.74 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Sylvie Laquerre, et al. 2009, EORTC International Conference. Abst B88.

Note: This product is for in vitro research use only