

FEN-1 Polyclonal Antibody

Catalog No: #40922

Package Size: #40922-1 50ul #40922-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FEN-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	FEN-1 Polyclonal Antibody detects endogenous levels of FEN-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human FEN-1.
Target Name	FEN-1
Other Names	FEN1; RAD2; Flap endonuclease 1; FEN-1; DNase IV; Flap structure-specific endonuclease 1; Maturation factor 1; MF1; hFEN-1
Accession No.	Swiss-Prot: P39748NCBI Gene ID: 2237
Uniprot	P39748
GeneID	2237;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using FEN-1 Polyclonal Antibody



Western Blot analysis of various cells using FEN-1 Polyclonal Antibody

Note: This product is for in vitro research use only