


Description

Product Name	L-701324
Brief Description	Inhibitors
Purification	99.50%
Target Name	NMDAR
Calculated MW	363.8
Formulation	C21H14ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC1=C(C(=O)NC2=C1C=CC(Cl)=C2)C1=CC(OC2=CC=CC=C2)=CC=C1
Pathways: Neuroscience
Receptor: NMDA receptor
Boiling pt: Melting pt: Solubility: DMSO: > 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wla? P, Poleszak E. Pharmacol Rep. 2011;63(5):1231-4.

Note: This product is for in vitro research use only