

Filamin 1 Polyclonal Antibody

Catalog No: #40927

Package Size: #40927-1 50ul #40927-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Filamin 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Filamin 1 Polyclonal Antibody detects endogenous levels of Filamin 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Filamin 1 around the non-phosphorylation site of S2152.
Target Name	Filamin 1
Other Names	FLNA; FLN; FLN1; Filamin-A; FLN-A; Actin-binding protein 280; ABP-280; Alpha-filamin; Endothelial actin-binding protein; Filamin-1; Non-muscle filamin
Accession No.	Swiss-Prot: P21333NCBI Gene ID: 2316
Uniprot	P21333
GeneID	2316;
SDS-PAGE MW	280kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Filamin 1 Polyclonal Antibody

Note: This product is for in vitro research use only