

Flt3 Polyclonal Antibody

Catalog No: #40931

Package Size: #40931-1 50ul #40931-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Flt3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	Flt3 Polyclonal Antibody detects endogenous levels of Flt3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Flt3 around the non-phosphorylation site of Y599.
Target Name	Flt3
Other Names	FLT3; CD135; FLK2; STK1; Receptor-type tyrosine-protein kinase FLT3; FL cytokine receptor; Fetal liver kinase-2; FLK-2; Fms-like tyrosine kinase 3; FLT-3; Stem cell tyrosine kinase 1; STK-1; CD antigen CD135
Accession No.	Swiss-Prot: P36888NCBI Gene ID: 2322
Uniprot	P36888
GeneID	2322;
SDS-PAGE MW	117kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HEPG2 cells using Flt3 Polyclonal Antibody

Note: This product is for in vitro research use only