

FoxO1A Polyclonal Antibody

Catalog No: #40936

Package Size: #40936-1 50ul #40936-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxO1A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	FoxO1A Polyclonal Antibody detects endogenous levels of FoxO1A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human FoxO1A around the non-phosphorylation site of S329.
Target Name	FoxO1A
Other Names	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma
Accession No.	Swiss-Prot: Q12778NCBI Gene ID: 2308
Uniprot	Q12778
GeneID	2308;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using FoxO1A Polyclonal Antibody



Western Blot analysis of HELA Jurkat cells using FoxO1A Polyclonal Antibody

Note: This product is for in vitro research use only