

FoxO3A Polyclonal Antibody

Catalog No: #40937

Package Size: #40937-1 50ul #40937-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxO3A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	FoxO3A Polyclonal Antibody detects endogenous levels of FoxO3A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human FoxO3A around the non-phosphorylation site of S253.
Target Name	FoxO3A
Other Names	FOXO3; FKHL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1
Accession No.	Swiss-Prot: O43524NCBI Gene ID: 2309
Uniprot	O43524
GeneID	2309;
SDS-PAGE MW	71kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using FoxO3A Polyclonal Antibody



Western Blot analysis of MCF7 cells using FoxO3A Polyclonal Antibody

Note: This product is for in vitro research use only