

Fra-1 Polyclonal Antibody

Catalog No: #40938



Package Size: #40938-1 50ul #40938-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Fra-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Fra-1 Polyclonal Antibody detects endogenous levels of Fra-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Fra-1.
Target Name	Fra-1
Other Names	FOSL1; FRA1; Fos-related antigen 1; FRA-1
Accession No.	Swiss-Prot: P15407NCBI Gene ID: 8061
Uniprot	P15407
GeneID	8061;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using Fra-1 Polyclonal Antibody

Note: This product is for in vitro research use only