

Frizzled-7 Polyclonal Antibody

Catalog No: #40939

Package Size: #40939-1 50ul #40939-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Frizzled-7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms
Specificity	Frizzled-7 Polyclonal Antibody detects endogenous levels of Frizzled-7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Frizzled-7.
Target Name	Frizzled-7
Other Names	FZD7; Frizzled-7; Fz-7; hFz7; FzE3
Accession No.	Swiss-Prot: O75084NCBI Gene ID: 8324
Uniprot	O75084
GeneID	8324;
SDS-PAGE MW	64kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HELA cells using Frizzled-7 Polyclonal Antibody



Western Blot analysis of HELA cells using Frizzled-7 Polyclonal Antibody

Note: This product is for in vitro research use only