

FRP-2 Polyclonal Antibody

Catalog No: #40940

Package Size: #40940-1 50ul #40940-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FRP-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	FRP-2 Polyclonal Antibody detects endogenous levels of FRP-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human FRP-2.
Target Name	FRP-2
Other Names	SFRP2; FRP2; SARP1; FKSG12; Secreted frizzled-related protein 2; FRP-2; sFRP-2; Secreted apoptosis-related protein 1; SARP-1
Accession No.	Swiss-Prot: Q96HF1NCBI Gene ID: 6423
Uniprot	Q96HF1
GeneID	6423;
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using FRP-2 Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN cells using FRP-2
Polyclonal Antibody

Note: This product is for in vitro research use only