

Fyn Polyclonal Antibody

Catalog No: #40943

Package Size: #40943-1 50ul #40943-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fyn Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Fyn Polyclonal Antibody detects endogenous levels of Fyn protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Fyn around the non-phosphorylation site of Y530.
Target Name	Fyn
Other Names	FYN; Tyrosine-protein kinase Fyn; Proto-oncogene Syn; Proto-oncogene c-Fyn; Src-like kinase; SLK; p59-Fyn
Accession No.	Swiss-Prot: P06241NCBI Gene ID: 2534
Uniprot	P06241
GeneID	2534;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Fyn Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Fyn Polyclonal Antibody

Note: This product is for in vitro research use only