

G2A Polyclonal Antibody

Catalog No: #40944

Package Size: #40944-1 50ul #40944-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	G2A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Mk
Specificity	G2A Polyclonal Antibody detects endogenous levels of G2A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human G2A.
Target Name	G2A
Other Names	GPR132; G2A; Probable G-protein coupled receptor 132; G2 accumulation protein
Accession No.	Swiss-Prot: Q9UNW8NCBI Gene ID: 29933
Uniprot	Q9UNW8
GeneID	29933;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using G2A Polyclonal Antibody



Western Blot analysis of various cells using G2A Polyclonal Antibody

Note: This product is for in vitro research use only