

Galectin-8 Polyclonal Antibody

Catalog No: #40946



Package Size: #40946-1 50ul #40946-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Galectin-8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Ms
Specificity	Galectin-8 Polyclonal Antibody detects endogenous levels of Galectin-8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Galectin-8.
Target Name	Galectin-8
Other Names	LGALS8; Galectin-8; Gal-8; Po66 carbohydrate-binding protein; Po66-CBP; Prostate carcinoma tumor antigen 1; PCTA-1
Accession No.	Swiss-Prot: O00214NCBI Gene ID: 3964
Uniprot	O00214
GeneID	3964;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of 3T3 cells using Galectin-8 Polyclonal Antibody

Note: This product is for in vitro research use only